



Alib
medical supplies limit

CATALOGUE

2025

MRI (1.5T / 0.35T / 0.5T)



NeuMR 1.5T

☒ Magnetic field strength	1.5T
☒ Max. gradient amplitude (mT/m)	35
☒ MRI Type	Open
☒ Gradient slew rate (mT/m/mcs)	136
☒ Type of magnet	Permanent
☒ Flat Panel Detector	17 × 17 inch



Superstar 0.35T

☒ Magnetic field strength	0.35T
☒ Max. gradient amplitude (mT/m)	26
☒ MRI Type	Open
☒ Gradient slew rate (mT/m/mcs)	67
☒ Type of magnet	Permanent

CT (64Slices)



NeuViz 64 In / En CT

☒ Focal spot sizes	0.6x1.2mm / 1.1x1.2mm
☒ Method of Cooling	Oil to air
☒ Powering Rate	50KW
☒ kV range	80,100,120,140KV
☒ mA range	30-420mA
☒ Detector rows	32rows-64slices

CT (16Slices)



NeuViz 16 Classic CT

☒ Focal spot sizes	0.4x0.7mm / 0.6x1.3mm
☒ Method of Cooling	Oil to air
☒ Powering Rate	28KW
☒ kV range	80,100,120,140KV
☒ mA range	10-230mA
☒ Detector rows	24 rows-16slices

Injector

CT Dual



CT Single



MR



Angio/DSA



Syringe