



Alib
medical supplies limit

CATALOGUE


2025

Digital X-ray

Fixed Digital X-ray



SONTU300-Mars



SONTU100-RAD



MTP50-B

Technical Configuration

Items / Model	SONTU300-Mars	SONTU100-RAD	MTP50-B
Output Power	50kW/50KHz	50kW/50KHz	50kW/60KHz
KV Range	40~150kV		
mA Range	100~630mA		100~650mA
Tube Focus	0.6mm/1.2mm		
Flat Panel Detector	17 × 17 inch	14 × 17 inch	17 × 17 inch

Ceiling Digital X-ray




SONTU-Polaris 50C

Technical Configuration

Items / Model	SONTU-Polaris 50C
Output Power	50kW/25KHz
KV Range	40~150kV
mA Range	100~630mA
Tube Focus	0.6mm/1.2mm
Flat Panel Detector	17 × 17 inch

U-ARM Digital X-ray



SONTU-P50E



DTP580B



Technical Configuration

Items / Model	SONTU-P50E	DTP580B
Output Power	50kW/50KHz	50kW/60KHz
KV Range	40~150kV	
mA Range	100~630mA	100~650mA
Tube Focus	0.6mm/1.2mm	
Flat Panel Detector	17 × 17 inch	

Mobile Digital X-ray





DP326

Technical Configuration

Items / Model	DP326
Output Power	32kW/50KHz
KV Range	40~150kV
mA Range	50~400mA
Tube Focus	0.6mm/1.2mm
Flat Panel Detector	14 × 17 inch

DRF Series





DTP576

Technical Configuration

Items / Model	DTP576
Voltage for Radiography	40kV-150kV
Voltage for Fluoroscopy	40kV-120kV
Current for Radiography	10mA-650mA
Current for Fluoroscopy	0.5mA-10mA
Tube Focus	0.6mm/1.2mm
SID	1000-1800mm
Flat Panel Detector	17 × 17 inch