



Alib
medical supplies limit

CATALOGUE

2025

Vehicle-mounted DR



Technical Configuration

Output Power	50kW
KV Range	40~150kV
mA Range	500mA
Tube Focus	0.6mm/1.2mm
Flat Panel Detector	17 × 17 inch
Tube Focus	0.6mm/1.2mm
SID	1000-1800mm
Flat Panel Detector	17 × 17 inch


Mobile C-arm System



PLX7000A/B/C

PLX7000/PLX118F Series
Technical Configuration

Items / Model	PLX7000A	PLX7000B	PLX7000C	PLX118F
Power/Frequency	6kW/60KHz	12kW/60KHz	16kW/60KHz	5kW/110KHz
Radiography Voltage	40kv~125kv			40kv-120kv
Radiography Current	120mA	160mA	200mA	100mA
Fluoroscopy Voltage	40kv-120kv			
Fluoroscopy Current	0.3mA~8mA	0.3mA~30mA		
Tube Focus	0.3mm/0.6mm			0.3/1.5mm
Image Intensifier	TOSHIBA (9")			Flat Panel



PLX112



PLX116

PLX112/PLX116 Series
Technical Configuration

Power/Frequency	5kW/40KHz
Radiography Voltage	40kv~120kv
Radiography Current	25-100mA
Fluoroscopy Voltage	40kv-120kv
Fluoroscopy Current	0.3mA~30mA
Tube Focus	0.6mm/1.8mm
Image Intensifier	TOSHIBA (9")