



**Alib**  
medical supplies limit

# CATALOGUE

2025



Digital Mammography



DM156  
Technical Configuration

| Items / Model       | DM156              |
|---------------------|--------------------|
| Output Power        | 5kW                |
| Frequency           | 40KHZ              |
| Output Voltage      | 20 - 35kV          |
| Cathodic Current    | 70-140mA / 20-30mA |
| Tube Focus          | 0.1mm/1.3mm        |
| SID                 | 650mm              |
| Flat Panel Detector | 10 × 12 inch       |


Analogue Mobile X-ray



PLX101 Series  
Technical Configuration

| Items / Model | PLX101        | PLX101C          | PLX101D |
|---------------|---------------|------------------|---------|
| Output Power  | 2.5kW/50KHz   | 5kW/50KHz        |         |
| KV Range      | 40~100kV      | 40-120kV         |         |
| mA Range      | 50mA          | 100mA            |         |
| Battery       | 220V, 50/60Hz | Built-in Battery |         |

Analogue Fixed X-ray



Technical Configuration

| Items / Model | YJF20R-A  | YJF32R-A | YJF50R-A  |
|---------------|-----------|----------|-----------|
| Output Power  | 20kW      | 32kW     | 50kW      |
| KV Range      | 40~125kV  | 40-125kV | 40-150kV  |
| mA Range      | 10-320mA  | 32-400mA | 10-630mA  |
| Focus         | 1.0/2.0mm |          | 0.6/1.2mm |