



**Alib**  
medical supplies limit

# CATALOGUE

2025

Stretcher / Examination Table

SKB039 (E)



Technical Parameters

|                |                |
|----------------|----------------|
| ☒ Full Size    | 1900×540×800mm |
| ☒ Lowest Size  | 1900×540×200mm |
| ☒ Back-rest    | 85°            |
| ☒ Net Weight   | 45kg           |
| ☒ Load Bearing | 159kg          |

X09



Technical Configuration

|                                        |      |
|----------------------------------------|------|
| ☒ Anti-slip Bed Feet Cover             | 4pcs |
| ☒ Mechanical Lever To Adjust Back Rest | 1set |

SKB040 (A)



Technical Parameters

|                |                |
|----------------|----------------|
| ☒ Full Size    | 1950×540×880mm |
| ☒ Lowest Size  | 1950×540×480mm |
| ☒ Net Weight   | 35kg           |
| ☒ Load Bearing | 159kg          |

X07a



Technical Configuration

|                            |      |
|----------------------------|------|
| ☒ Anti-slip Bed Feet Cover | 4pcs |
|----------------------------|------|

Gynecological Examination Table

A045-1



Technical Configuration

|                            |      |
|----------------------------|------|
| ☒ Lithotomy Rod            | 1set |
| ☒ Hand Bar                 | 1set |
| ☒ Anti-slip bed feet cover | 4pcs |
| ☒ Mechanical lever         | 2pcs |
| ☒ Stainless steel basin    | 1pcs |

Gynecological Examination Table

A99-8



Technical Configuration

|                            |       |
|----------------------------|-------|
| ☒ Hand Bar                 | 2 pcs |
| ☒ Leg Support With Cushion | 2 pcs |
| ☒ Foot Pedal               | 1pcs  |
| ☒ Stainless Steel Basin    | 1pcs  |